




# Monitoring Electrodermal Activity in RapidABC Sessions

Ivan Riobo<sup>1</sup>, Javier Hernandez<sup>2</sup>, Gregory D. Abowd<sup>1</sup>, Rosalind W. Picard<sup>2</sup>  
{ivan.riobo,abowd}@gatech.edu – Georgia Tech<sup>1</sup>, {javierhr, picard}@media.mit.edu – MIT Media Lab<sup>2</sup>

Emotional regulation and co-regulation are key components of engagement. Measuring EDA during RapidABC sessions may provide insightful information about social interactions with children.

## Electrodermal Activity (EDA)



### Good Indicator

- Arousal
- Cognitive Load
- Emotions

### Wearable Sensors

- Wireless
- Comfortable

### Very Sensitive

- Temperature
- Physical Activity
- Pressure

### Measurement Locations

- Wrists
- Ankles

(also known as Galvanic Skin Response)

## Data Collection Characteristics

- Sampling Rate** 32 Hz
- Sensors used per child** 2 (1 per wrist)
- Electrodes per sensor** 2
- Electrode Gel to enhance conductivity**
- Main Challenges**
  - Obtaining EDA baselines
  - Individual differences
  - Motion and contact artifacts
  - Synchronization between sensors
  - Synchronization with other modalities

## EDA Responses in RapidABC

### Positive Arousal

### Negative Arousal

Right and Left EDA

